



THE GLOBAL LEADER IN GEL TECHNOLOGIES®

Product Catalog

orthopedics

prosthetics

skin care



orthopedics

prosthetics

skin care

Our story...

Founded in 1989, Silipos has established itself as the global leader in gel technologies across a wide range of industries including healthcare, consumer goods and health and beauty aids. Our products are manufactured in the United States and sold to over 1,000 customers in 60 countries worldwide.

Silipos is committed to expanding and improving our range of products to provide our customers with new and innovative gel solutions. We pride ourselves on product excellence and quality, and are pleased to announce that we are an ISO 9001:2000 and ISO 13485:2003 certified company. Additionally, many of our products are CE marked for distribution in worldwide markets.

Within the healthcare industry, Silipos is the most trusted source for medical solutions among healthcare professionals. We have built our reputation on innovation, quality and outstanding customer service.

Some of the medical specialties to which we provide service are:

- Podiatry
- Orthopedics
- Prosthetics
- Physical Therapy
- Occupational Therapy
- Sports Medicine
- Scar Management
- Plastic Surgery
- Dermatology
- Cosmetology
- Rheumatology
- Home Healthcare

Our ongoing mission is to provide exemplary products and services that improve the life quality of our customers and community. Silipos' mineral oil gel products offer unique benefits unavailable in the competitive market.

Silipos' sales and marketing office is located in New York City, New York and our manufacturing facility is located in Niagara Falls, New York.



orthopedics 6

Digital Care	7
Bunion Protection	14
Metatarsal	16
Sleeves & Wraps	20
Insoles & Socks	24
Gel Pads & Sheeting	30
TheraStep®	34

prosthetics 41

Sheaths	42
Liners	44
Sleeves	52
Padding & Accessories	54

skin care 56

Gel Gloves & Socks	57
Scar Management	58



Our gels...

Silipos gels cushion and protect specific parts of the body against impact, pressure, and shear forces, while simultaneously releasing moisture-rich mineral oil that gently hydrates, softens, and conditions the skin without the use of additional creams or lotions.

Visco-elastic properties enable our gels to move synergistically with the skin which prevents abrasion and discomfort. Our user-friendly products easily adapt to a wide range of body types and builds to ensure high patient compliance.

All of our gels are clear, non-toxic, and hypoallergenic. Technologically advanced manufacturing capabilities allow Silipos to create a vast assortment of gels for a variety of uses. Additives such as vitamins and minerals, antioxidants, moisturizers, and fragrances can be applied to our gels to address specific concerns across numerous markets. Silipos gels are washable and reusable, and with proper care will continue to release skin-beneficial mineral oils throughout the life of the product.



Orthopedics

Silipos gels easily conform to the body ensuring a comfortable fit that cushions and protects. Made with mineral oil gels, our orthopedic and foot care products provide cushioning protection with the added benefit of skin moisturizing and conditioning when the gel is applied directly to the skin.

Our orthopedic products offer easy to use solutions for a wide variety of ailments. From common foot discomforts to post traumatic injuries, Silipos has a product that can help improve the quality of life for your patients.



Digital Pads

- Protects toes and fingers from friction and abrasion
- Helps to alleviate footwear pressure on toes
- Use where corns, calluses, or hammertoes cause discomfort
- Packaged in a polybag



Item #	Size	Quantity
10055	Mini	12/pkg
10065	S/M	12/pkg
10075	L/XL	12/pkg



Digital Pads on a Strip

- Mineral oil gel pads on a continuous strip
- Elastic tubing is easy to cut to preferred size
- Once cut, product is turned inside out so gel can be applied to skin
- Provides a more cost effective alternative to Digital Pads
- Packaged in a polybag

Item #	Size	Quantity
1105	Mini	17 pads per strip, 2 strips/pkg
1106	S/M	16 pads per strip, 2 strips/pkg
1107	L/XL	14 pads per strip, 2 strips/pkg



DigiStrip


- Half coated with mineral oil gel strip
- Elastic tubing is easy to cut to preferred size
- Once cut, product is turned inside out so gel can be applied to skin
- Provides a more cost effective alternative to our Corn Pads
- Approximately 30 pads per strip (when cut to 3cm lengths)
- Packaged in a polybag

Item #	Size	Quantity
1100	Narrow, 24" (60cm) long	2 strips/pkg
1101	Medium, 24" (60cm) long	2 strips/pkg
1102	Wide, 24" (60cm) long	2 strips/pkg



Corn Pads

- Helps to relieve painful corns and calluses
- Each pad is 3cm in length and contains 2mm thick mineral oil gel pad
- Gel pad allows for custom placement
- Packaged in a polybag




Item #	Size	Quantity
90055	Narrow	100/pkg
90065	Medium	100/pkg
90075	Wide	100/pkg



Gel Tubing

- Fully lined with mineral oil gel
- Provides all around digital protection
- Not as bulky as conventional foam tubing
- Durable and will not bottom out as foam materials do
- Can aid in reducing scar tissue
- Cut to desired size
- Packaged in a polybag



Ribbed Knit Fabric



Item #	Size	Quantity
10865	Narrow, 3/4" (2cm) diameter, 6" (15cm) length	2/pkg
10875	Wide, 1" (2.5cm) diameter, 6" (15cm) length	2/pkg

Mesh Fabric (for greater compression)



Item #	Size	Quantity
10855	Narrow, 3/4" (2cm) diameter, 6" (15cm) length	3/pkg
10845	Wide, 3" (7.6cm) diameter, 10" (25.4cm) length	1/pkg



Gel Toe Crest

- Helps to relieve pressure and pain caused by hammer, claw, and mallet toes
- Silicone gel designed to absorb shock and provide even pressure distribution
- Reinforced toe loop to ensure proper fit
- Will not bottom out like foam products
- Available in polybag or box packaging

Polybag Packaging

Item #	Size	Quantity
10415	Small, Left	3/pkg
10425	Small, Right	3/pkg
10435	Large, Left	3/pkg
10445	Large, Right	3/pkg

Box Packaging

Item #	Size	Quantity
10515	Small, Left	1/pkg
10525	Small, Right	1/pkg
10535	Large, Left	1/pkg
10545	Large, Right	1/pkg



Digital Caps

- Fully lined with mineral oil gel for all-around digital protection
- Helps protect against pain caused by pressure and friction
- Can be used to help treat nail or skin problems due to dryness
- Helps reduce discomfort caused by corns, blisters, and ingrown nails
- Fits toes or fingers (average cap length is 2.5"/6.5cm) and can be cut to size
- Packaged in a polybag



Ribbed Knit Fabric



Item #	Size	Quantity
10195	Mini	6/pkg
10205	S/M	6/pkg
10206	L/XL	6/pkg
10207	XXL	6/pkg

Mesh Fabric (for greater compression)



Item #	Size	Quantity
10235	S/M	8/pkg



Digital Cap with Toe Spreader

- Traditional Digital Cap with added benefit of built-in Gel Toe Spreader
- Fully lined with mineral oil gel for all-around digital protection
- Helps protect against pain caused by pressure and friction
- Comfortable and simple solution for patients suffering from multiple ailments
- Packaged in a polybag



Item #	Size	Quantity
1596	L/XL	2/pkg



All Gel Digital Caps

- Helps protect against pain caused by pressure and friction
- Can be used to help treat nail or skin problems due to dryness
- Helps reduce discomfort caused by corns, blisters, and ingrown nails
- Wear with sandals or open-toe shoes
- Packaged in a polybag




Item #	Size	Quantity
10495	S/M	2/pkg
10505	L/XL	2/pkg



Cushion Digital Cap

- Soft breathable foam cap provides lightweight comfort and protection
- Helps reduce discomfort caused by corns, blisters, and ingrown nails
- Packaged in a polybag



Item #	Size	Quantity
1589	S/M	4/pkg
1590	L/XL	4/pkg



Gel Toe Separators

- Anatomically designed to fit comfortably between toes
- Relieves pressure and friction between toes
- Helps correct overlapping toes
- Packaged in a polybag

Item #	Size	Quantity
11605	Small	15/bag
11615	Large	15/bag



Cushion Toe Separators

- Relieves pressure and friction between toes
- Helps correct overlapping toes
- Soft breathable foam provides lightweight comfort and protection
- Packaged in a polybag

Item #	Size	Quantity
1584	Small	6/pkg
1585	Medium	6/pkg
1586	Large	6/pkg



Gel Toe Spreaders

- Conforms and comfortably fits between toes
- Larger in size compared to the Gel Toe Separators
- Used as an inter-digital wedge for bunion deformity, overlapping toes, or toe drift
- Packaged in a polybag


Item #	Size	Quantity
11505	Small	4/bag
11515	Medium	4/bag
11525	Large	4/bag





Cushion Toe Spreaders

- Conforms and comfortably fits between toes
- Helps relieve pain caused by bunions
- Soft breathable foam provides lightweight comfort and protection
- Packaged in a polybag




Item #	Size	Quantity
1587	Medium	4/pkg
1588	Large	4/pkg



Cushion Toe Comb

- Separates, cushions, and protects toes
- Helps relieve discomfort caused by crooked or overlapping toes
- Soft breathable foam provides lightweight comfort and protection
- Packaged in a polybag



Item #	Size	Quantity
1591	One Size	6/pkg



Digital Care Kit

- A convenient and economical way to store, display, and dispense our most popular digital care products
- Special compartments for each product allows you to easily select the best solution for the condition you are treating

Kit Includes:



Digital Caps

- (6) S/M

Digital Pads

- (4) Mini
- (4) S/M
- (4) L/XL

Gel Toe Separators

- (6) Small
- (6) Large

Gel Toe Spreaders

- (4) Medium
- (4) Large

Mesh Tubing

- (2) 3/4" x 6"

Item #	Quantity
30365	1 Kit



All Gel Care Kit

- Offers an assortment of our most popular gel foot care products in a compact kit for easy in-office dispensing
- An excellent compliment to the original Digital Care Kit

Kit Includes:



All Gel Digital Caps

- (4) S/M
- (4) L/XL

Gel Toe Separators

- (6) Small
- (6) Large

Gel Metatarsal Pad w/ Toe Spreader

- (Pair) S/M
- (Pair) L/XL

Gel Toe Spreaders


- (4) Medium
- (4) Large

Item #	Quantity
31744	1 Kit



Bunion Care Gel Sleeve

- Anatomically designed to comfortably conform to the shape of the bunion
- Mineral oil gel pad cushions, protects, and reduces pressure on affected area
- Helps to reduce scarring after surgery
- Packaged in a polybag




Item #	Size	Quantity
10325	S/M	1/pkg
10335	L/XL	1/pkg



Gel Bunion

- Anatomically designed to comfortably conform to the shape of the bunion
- Mineral oil gel pad cushions, protects, and reduces pressure on affected area
- Helps to reduce scarring after surgery
- Packaged in a plastic clamshell




Item #	Size	Quantity
10085	One Size	5/pkg



Slim Gel-Fit™ Bunion Sleeve

- Elastic bunion sleeve with the gel pad incorporated inside two layers of a comfortable and conforming fabric
- The covered gel pad shields the painful and sensitive bony protrusion of the forefoot
- Product helps relieve shoe pressure, friction, and shear forces on the bunion
- May also be worn on the hand to protect the thumb (thenar eminence) from discomfort during gripping activities
- Product's design is thin and compact making it perfect for everyday use and fit with most dress, casual, or athletic shoes and underneath hosiery
- Packaged in a polybag




Item #	Size	Quantity
1110	Small	1/pkg
1111	Medium	1/pkg



Gel Foot Cover

- Thin layer of gel covered with fabric offers complete toe protection
- Extends to provide protection and cushioning to the metatarsal heads
- Designed to be worn comfortably with everyday footwear and underneath hosiery
- Ideal for activities that place stress on toes and metatarsals
- Packaged in a polybag




Item #	Size	Quantity
10275	Small, 3.75" (9.5cm) x 4.9" (12.5cm)	1 pair/pkg
10285	Medium, 4" (10cm) x 5.3" (13.5cm)	1 pair/pkg
10295	Large, 4.3" (11cm) x 5.9" (15cm)	1 pair/pkg



Universal Gel Strap

- Comfortable elastic sleeve with anatomically designed gel pad to cover all five metatarsal heads
- Helps reduce pain and discomfort in the forefoot
- Serves as a soft tissue supplement for fat pad atrophy
- Multi-functional use; cut stitching to fit to instep, bunions, elbow, hand, or wrist
- Available with uncovered (exposed) gel pad or covered gel pad
- Packaged in a polybag



Gel Uncovered

Mineral oil gel pad helps soften and reduce callused skin

Item #	Size	Quantity
10155	S/M	1 pair/pkg
10165	L/XL	1 pair/pkg

Gel Covered


Extended product life when callusing is not an issue

Item #	Size	Quantity
10175	S/M	1 pair/pkg
10185	L/XL	1 pair/pkg



Slim Gel-Fit™ Metatarsal Pad

- Elastic gel metatarsal sleeve with the gel pad incorporated inside two layers of a comfortable and conforming fabric
- The covered gel pad eases pain under the ball of the foot caused by metatarsalgia, calluses, rheumatism, neuromas, or from long hours on your feet
- Product's design is thin and compact making it perfect for everyday use and fit with most dress, casual, or athletic shoes and underneath hosiery
- Great for use with high heels that put extra pressure on metatarsal heads
- Packaged in a polybag



Item #	Size	Quantity
1113	Small	1 pair/pkg
1114	Medium	1 pair/pkg



Ball of Foot Gel Cushion

- Thin fabric-covered gel pad (2mm) to protect against abrasion, friction, and pressure
- Excellent for relieving the burning sensation under the metatarsal heads
- May be used as a soft tissue supplement for metatarsalgia
- Ideal for calluses and small neuromas
- Thin enough to fit in any shoe type
- Packaged in a polybag

Item #	Size	Quantity
10025	One Size	2/pkg



Gel Metatarsal Pad

- Absorbs shock, impact, and provides even distribution of pressure in the forefoot area
- Anatomically designed for comfort and support
- Helps to reduce the burning sensation in the metatarsal heads
- Ideal for conditions such as fat pad atrophy, metatarsalgia, and Morton's neuroma
- Silicone construction with reinforced toe loop to ensure proper fit
- Packaged in a polybag

Item #	Size	Quantity
10465	One Size	1 pair/pkg



Gel Metatarsal Pad with Toe Spreader

- Anatomically designed for comfort and protection of metatarsal heads
- Helps reduce pain and burning sensations under metatarsal heads
- Ideal for fat pad atrophy, metatarsalgia, and Morton's neuroma
- Silicone construction with built-in toe spreader helps reduce discomfort between toes while securing product in place
- Ideal for activities that place stress on toes and metatarsals
- Packaged in a box

Item #	Size	Quantity
4020	S/M	1 pair/pkg
4021	L/XL	1 pair/pkg



Gel Ball of Foot Pad

- Absorbs shock and pressure on the ball of foot
- Gently massages the bottom of the foot
- Helps prevent calluses by reducing friction and abrasion
- Packaged in a polybag

Item #	Size	Quantity
1582	One Size	1 pair/pkg



Cushion Ball of Foot Pad


- Absorbs shock and pressure on the ball of foot
- Gently massages the bottom of the foot
- Helps prevent calluses by reducing friction and abrasion
- Packaged in a polybag

Item #	Size	Quantity
1583	One Size	1 pair/pkg



Achilles Heel Pad

- Soft conforming gel pad protects the posterior of the heel from pressure, impact, and friction
- Provides gentle compression and dissipation of pressure around inflamed area
- Effective in protecting against discomfort caused by ski boots, skates, cleats, other athletic footwear, and shoes or boots
- Reduces pain caused by footwear pressure on irregular heel bumps
- Packaged in a polybag




Item #	Size	Quantity
10385	S/M	1/pkg
10395	L/XL	1/pkg



Malleolar Sleeve

- Designed with two discs of gel molded on both the medial and lateral malleolae
- Cushions and protects the ankle against impact, pressure, and abrasion
- Provides gentle compression to the ankle
- Effective in protecting against discomfort caused by ski boots, skates, cleats and other athletic footwear
- Packaged in a polybag




Item #	Size	Quantity
15005	One Size	1/pkg



Soft Skin Heel Sleeve

- Thin gel pad anatomically molded to fit the heel
- Designed to moisturize and soften dry, cracked skin and help to reduce heel fissures
- Beneficial for dermatitis and eczema
- Can be used in conjunction with skin moisturizing creams and lotions
- Recommended as a non-ambulatory therapy product, not to be worn with shoes
- Packaged in a polybag





Item #	Size	Quantity
15265	One Size	1/pkg



Heel/Elbow Slip Over

- Versatile mineral oil gel pad can be positioned almost anywhere around the foot and heel
- Aids in the healing and skin conditioning of dry cracked heels and elbows
- Effective in helping to prevent heel and elbow ulcerations
- Packaged in a polybag


Item #	Size	Quantity
15225	S/M	1/pkg
15235	L/XL	1/pkg



Boot Bumper

- Helps protect the foot and ankle from abrasions without the added bulk
- Effective in protecting against discomfort caused by ski boots, skates, cleats and other athletic footwear
- Provides light compression to the foot and ankle
- Packaged in a plastic clamshell


Item #	Size	Quantity
1770	S/M	1/pkg
1772	L/XL	1/pkg



Skate Bite Protector

- Anatomically designed gel pad to protect the dorsal aspect of the foot from pressure, friction, and pinching from shoelaces
- Effective in protecting against discomfort caused by ski boots, skates, cleats and other athletic footwear
- Provides light compression to the foot and ankle
- Packaged in a polybag

Item #	Size	Quantity
92573	One Size	1/pkg






Arch Support

- Elastic sleeve provides gentle compression and comfort to arch of foot
- Comfortably lifts and supports the plantar fascia
- Helps to relieve heel and arch pain
- Packaged in a polybag

Item #	Size	Quantity
1595	One Size	2/pkg



Half Coated Body Sleaving

- New and improved fabric for better fit and flexion
- Knit sleaving with moisturizing mineral oil gel coating on one side of the sleeve only
- Available in 10 foot rolls allowing sleeve to be dispensed and cut to desired length
- Excellent for use on upper and lower extremity allowing for many universal uses
- Moisturization and mild compression attributes make this product excellent for: skin grafts, burns, scar treatment, and Lymphedema applications
- Packaged in a polybag

Item #	Size	Quantity
10925	3" (8.3cm) x 10' (3.1m)	1 roll/pkg
10935	4" (10cm) x 10' (3.1m)	1 roll/pkg
10945	5" (13cm) x 10' (3.1m)	1 roll/pkg



Carpal Gel Sleeve

- Anatomically designed gel pad molded to elastic sleeve contours to shape of palm
- Helps prevent and reduce shock, impact, and vibration
- Reduces discomfort caused by hypersensitivity
- For placement directly over carpal tunnel post-operative scar area
- Mineral oil gel pad helps prevent and reduce scar tissue
- Can be worn under gloves and/or splints
- Packaged in a polybag

Item #	Size	Quantity
14115	Small, Left	1/pkg
14125	Small, Right	1/pkg

Item #	Size	Quantity
14135	Large, Left	1/pkg
14145	Large, Right	1/pkg



Gel-E-Roll

- Elastic bandage with one side coated in a thin layer of mineral oil gel
- Conforms to difficult or irregular surfaces of the body and will not slide
- Gel moves with the skin and does not promote skin irritation
- Packaged in a polybag
- May be used:

- to reduce scar tissue
- for cold therapy
- to provide compression
- as a residual limb shrinker
- to aid in muscle strains or arthritic pain
- to moisturize areas of dry skin
- for burn care
- to help control edema

Item #	Size	Quantity
10355	4" (10cm) x 48" (120cm)	1/pkg

WonderZorb®

The #1 off-the-shelf foot orthotic prescribed by healthcare professionals for treating foot pain

- Dual density orthotics provide maximum support and shock absorption.
- Clear silicone construction with a softer density colored insert provides additional relief to sensitive areas.



WonderSpur

- Aids in calcaneal spurs, plantar fasciitis, heel spurs, and knee & back pain
- Packaged in a box

Item #	Size	Mens	Womens	Quantity
4500	Small	2-4	4-5	1 pair/pkg
4501	Medium	5-7	6-9	1 pair/pkg
4502	Large	8-10	10-12	1 pair/pkg
4503	X-Large	11+	12+	1 pair/pkg



WonderCup

- Relieves pain throughout heel
- Aids in bursitis, achilles tendonitis, pressure ulcers, and plantar fasciitis
- Packaged in a box

Item #	Size	Mens	Womens	Quantity
4504	Small	2-4	4-5	1 pair/pkg
4505	Medium	5-7	6-9	1 pair/pkg
4506	Large	8-10	10-12	1 pair/pkg
4507	X-Large	11+	12+	1 pair/pkg



WonderSport

- Simultaneously cushions and supports heel, and arch
- Aids in fat pad atrophy and with diabetic feet
- Packaged in a box

Item #	Mens	Quantity
4508	5-7	1 pair/pkg
4509	8-9	1 pair/pkg
4510	10-11	1 pair/pkg
4511	12-13	1 pair/pkg
4512	14+	1 pair/pkg

Item #	Womens	Quantity
4513	4-5	1 pair/pkg
4514	6-7	1 pair/pkg
4515	8-9	1 pair/pkg
4516	10+	1 pair/pkg



WonderSole

- Full length insole designed to provide support to metatarsals, arch, and heel
- Helps relieve symptoms related to metatarsalgia, plantar fasciitis, fat pad atrophy, and arthritis
- Packaged in a box

Item #	Size	Mens	Womens	Quantity
4530	Small	7	5-7	1 pair/pkg
4531	Medium	8-9	8-9	1 pair/pkg
4532	Large	10	10-11	1 pair/pkg
4533	X-Large	11-13	—	1 pair/pkg



3/4 Length Insole

- Anatomically designed to support longitudinal and transverse arches
- Deep heel cup ensures superior comfort and control
- Soft dual-density silicone provides shock absorption and pressure relief to sensitive areas
- Packaged in a box

Item #	Size	Mens	Womens	Quantity
4010	Small	—	4.5-6	1 pair/pkg
4011	Medium	5-7.5	6.5-9.5	1 pair/pkg
4012	Large	8-10	10-11.5	1 pair/pkg



Full Length Insole

- Multi-density silicone absorbs shock under high pressure areas while supporting the rest of the foot
- Ideal as a replacement for shoes with removable insoles
- Relief for heel, arch, and forefoot problems
- Multiple vents enhance air circulation
- Packaged in a box

Item #	Size	Mens	Womens	Quantity
4013	Small	5-6.5	6.5-8.5	1 pair/pkg
4014	Medium	7-9.5	9-11	1 pair/pkg
4015	Large	10-12	11.5-12	1 pair/pkg
4016	X-Large	12.5	—	1 pair/pkg



Full Length Insole with Toe Crest

- Thin multi-density silicone design ensures full-foot comfort and protection and is ideal as a replacement for shoes with removable insoles
- Has soft full built-in toe crest for added functionality and metatarsal raise for relief of forefoot discomfort
- Designed to reduce pressure while supporting mid and forefoot areas
- Packaged in a box

Item #	Size	Mens	Womens	Quantity
4017	Small	—	5-6	1 pair/pkg
4018	Medium	5-6	6.5-7.5	1 pair/pkg
4019	Large	6.5-7.5	8-9	1 pair/pkg
4033	X-Large	8-9	10-10.5	1 pair/pkg
4034	XX-Large	9.5-10	11-11.5	1 pair/pkg



Low Profile Fabric-Covered Insole

- For general shock absorption, pressure distribution, and added comfort
- Soft, durable, and comfortable thin fabric-covered insole allows foot to slide easily into shoe
- May be used as a lining under other insoles for additional shock absorption and protection
- Packaged in a box

Item #	Size	Mens	Womens	Quantity
4030	Small	—	4.5-7	1 pair/pkg
4031	Medium	5.5-9.5	7.5-11	1 pair/pkg
4032	Large	10-12.5	11.5-12	1 pair/pkg



Heel Spur with Removeable Plug

- Shock absorbing silicone heel cushion with removable plug
- Wear with plug as a comfortable heel pad, wear without plug for heel relief and protection around the heel spur
- Fits all shoe styles and types, will not slide in shoe
- Packaged in a box

Item #	Size	Mens	Womens	Quantity
4001	Small	—	5-6.5	1 pair/pkg
4002	Medium	5-7.5	7-9.5	1 pair/pkg
4003	Large	8-11.5	10-12	1 pair/pkg



Low Profile Heel Cushion

- Soft silicone heel cushion absorbs shock and evenly distributes weight
- Low profile design fits comfortably in all shoe styles and types, will not slide in shoe
- Packaged in a box

Item #	Size	Mens	Womens	Quantity
4004	Small	—	5-6.5	1 pair/pkg
4005	Medium	5-7.5	7-9.5	1 pair/pkg
4006	Large	8-11.5	10-12	1 pair/pkg



Cushion Heel Pad

- Helps relieve heel pain
- Self adhesive backing for proper placement
- Reduces discomfort caused by plantar fasciitis
- Packaged in a polybag

Item #	Size	Quantity
1592	Small	1 pair/pkg
1593	Medium	1 pair/pkg
1594	Large	1 pair/pkg



Arthritic/Diabetic Gel Sock

- Single layer sock with 1/16" (2mm) of gel throughout the sole
- Helps prevent foot ulcerations
- Non-restrictive cuff
- Can be used with or without diabetic shoes
- Provides optimum protection for the neuropathic foot
- Helps reduce friction, abrasion, and shear forces while preventing callus formation
- Serves as a soft tissue supplement for fat pad atrophy
- Packaged in a box

☐ **White**

Item #	Size	Quantity
1701	US Size 7-9 (UK Size 37-40) Shoe: Womens 6 or smaller	1 pair/pkg
1702	US Size 9-11 (UK Size 41-43) Shoe: Womens 6 and up / Mens 9 or smaller	1 pair/pkg
1703	US Size 11-13 (UK Size 44-46) Shoe: Mens 9 and up	1 pair/pkg

■ **Black**

Item #	Size	Quantity
1706	US Size 7-9 (UK Size 37-40) Shoe: Womens 6 or smaller	1 pair/pkg
1707	US Size 9-11 (UK Size 41-43) Shoe: Womens 6 and up / Mens 9 or smaller	1 pair/pkg
1708	US Size 11-13 (UK Size 44-46) Shoe: Mens 9 and up	1 pair/pkg



SoftSock

- Double layer sock with 1/4" (6mm) of gel throughout the sole
- Excellent for any condition that requires a soft, friction absorbing foot bed
- Indicated for neuropathic & rheumatoid arthritis sufferers
- Fluid like properties of the gel help massage and comfort the foot while walking
- Effective in reducing shock and impact on the heel
- Serves as a soft tissue supplement for fat pad atrophy
- Requires extra depth footwear
- Packaged in a polybag

☐ **White**

Item #	Size	Quantity
10605	Small, US Size 7-9 (UK Size 37-40)	1 pair/pkg
10615	Medium, US Size 9-11 (UK Size 41-43)	1 pair/pkg
10625	Large, US Size 10-13 (UK Size 43-46)	1 pair/pkg
10635	X-Large, US Size 13-16 (UK Size 46+)	1 pair/pkg



Gel Dots

- Self-adhesive backing allows for use with most materials (not skin)
- Gel dot is 1/8" (3mm) thick and 1" (3cm) in diameter
- Packaged in a polybag
- May be used:
 - to reduce abrasion, friction, and pressure points in splints, casts, braces, or footwear
 - can be fabricated with insoles and placed under individual metatarsal heads
 - to protect the ulnar styloid
 - as individual cushions for parts of the hand such as the knuckles in casts or splints

Item #	Quantity
10255	15 dots/sheet



Gel Squares

- Versatile solution for discomfort related to pressure and friction
- Gel square is 1/8" (3mm) thick and 4" (10cm) x 4" (10cm) in size
- Can be trimmed to customized size with scissors
- Mineral oil gel moisturizes and conditions skin
- Can be used as an interface between the skin and a brace
- Packaged in a polybag

Item #	Size	Quantity
15505	With Self-Adhesive Backing	2 pieces/pkg
15515	Without Backing	2 pieces/pkg



Body Discs

- Multi-purpose discs of mineral oil gel
- Apply to parts of the body suffering from discomfort caused by pressure and friction
- Can be trimmed to customized size with scissors
- Mineral oil gel moisturizes and conditions skin
- Packaged in a polybag

Item #	Size	Quantity
10305	2.5" (6.3cm) Diameter	2 pieces/pkg
10315	4" (10cm) Diameter	2 pieces/pkg



Carpal Gel Pads

- Anatomically designed gel pad (same as in Carpal Gel Sleeve) with self-adhesive backing
- Can be incorporated into custom made or pre-fabricated hand and wrist braces
- Packaged in a polybag

Item #	Size	Quantity
14205	2" (5cm) x 6" (15.24cm)	3/pkg



Pressure Relief Padding

- 2mm thick gel pad with self-adhesive backing
- Used to reduce shock, impact, and vibration
- Can be used as a replacement for moleskin or felt products
- Excellent for use with custom ankle foot orthosis devices
- Perfect for cushion protection in hand splints
- Packaged in a cardboard tube

Item #	Size	Quantity
11205	4" (10cm) x 36" (91cm)	1/pkg
11215	8" (20cm) x 36" (91cm)	1/pkg



Soft Shear Gel Sheeting

- 1/8" (3mm) thick gel sheet available with mineral oil gel exposed on one side, or with both sides lined with fabric
- Fabric top cover provides excellent stretch properties
- Relieves pressure and shear forces
- Fabric side can be glued to braces or other materials
- Shrinkwrap packaged with backerboard

Fabric One Side

Item #	Size	Quantity
15605	8" (20cm) x 8" (20cm)	1/pkg
15705	12" (30.5cm) x 16" (40.5cm)	1/pkg
15805	16" (40.5cm) x 20" (51cm)	1/pkg

Fabric Two Sides

Item #	Size	Quantity
15615	8" (20cm) x 8" (20cm)	1/pkg
15715	12" (30.5cm) x 16" (40.5cm)	1/pkg
15815	16" (40.5cm) x 20" (51cm)	1/pkg



WonderFlex

- Silicone sheeting available in multiple thicknesses and in lined or unlined versions
- Wonderflex sheeting is 12" (30.5cm) x 16" (40.5cm) in size
- Shrinkwrap packaged with backerboard

Unlined

Item #	Gel Thickness	Quantity
4722	2mm	1/pkg
4723	3mm	1/pkg

- *Recommended for use with orthotic devices*

Lined

Item #	Gel Thickness	Quantity
4733	3mm	1/pkg
4735	5mm	1/pkg

- *Soft durometer cushioning for top cover of orthotics*
- *Ideal for diabetic, arthritic, and sensitive feet*
- *Resists compression*



Silipos is pleased to introduce the next step in foot therapy...

Cushion, protect and relieve foot pain and pressure with every step. TheraStep®'s comforting foot care products turn uncomfortable walks into soothing therapy sessions for troubled feet.

Each product relieves the pain associated with bunions, corns, blisters and calluses comfortably and effectively, as mineral oil conditioners moisturize and soften skin.

New cushion foam products provide a lightweight, breathable alternative.

The TheraStep® consumer packaged line offers the most comprehensive program for professional in-office dispensing. Increase revenue for your practice while providing your patients with a proven solution for their ailments.





Gel Toe Cap

- Cushioning protection for toes or fingers
- Helps reduce discomfort caused by corns, blisters, and ingrown nails
- Mineral oil gel soothes and moisturizes the skin

Item #	Size	Quantity
6203	Small	1/pkg
6204	Large	1/pkg



All-Gel Toe Cap

- Cushioning protection for toes or fingers
- Wear with sandals or open-toe shoes
- Helps reduce discomfort caused by corns, blisters, and ingrown nails

Item #	Size	Quantity
6201	Small	1/pkg
6202	Large	1/pkg



Cushion Toe Cap

- Cushioning protection for toes or fingers
- Helps reduce discomfort caused by corns, blisters, and ingrown nails
- Soft breathable foam provides lightweight comfort

Item #	Size	Quantity
6218	Small	1/pkg
6219	Large	1/pkg





Gel Corn Pads

- Helps to relieve painful corns
- Gel pad allows for custom placement
- Mineral oil gel soothes and moisturizes the skin

Item #	Size	Quantity
6205	One Size	2/pkg



Gel Toe Tubes

- Fully lined with gel for all-around protection
- Cushions and protects toes against pressure and friction
- Mineral oil gel soothes and moisturizes the skin
- Can be cut to size

Item #	Size	Quantity
6208	One Size	2/pkg



Gel Toe Separators

- Relieves pressure and friction between toes
- Helps correct overlapping toes
- Mineral oil gel soothes and moisturizes the skin

Item #	Size	Quantity
6206	2 Small, 2 Large	4/pkg





Cushion Toe Separators

- Relieves pressure and friction between toes
- Helps correct overlapping toes
- Soft breathable foam provides lightweight comfort

Item #	Size	Quantity
6216	2 Small, 2 Medium, 2 Large	6/pkg



Gel Toe Spreaders

- Helps relieve pain caused by bunions
- Helps correct overlapping toes
- Mineral oil gel soothes and moisturizes the skin

Item #	Size	Quantity
6207	Medium	2/pkg



Cushion Toe Spreaders

- Helps relieve pain caused by bunions
- Helps correct overlapping toes
- Soft breathable foam provides lightweight comfort

Item #	Size	Quantity
6217	1 Medium, 1 Large	2/pkg





Cushion Toe Comb

- Separates and cushions toes
- Helps relieve discomfort caused by crooked or overlapping toes
- Soft breathable foam provides lightweight comfort

Item #	Size	Quantity
6220	One Size	2/pkg



Gel Toe Support

- Helps relieve toe pain caused by hammer, claw, and mallet toes
- Reduces toe and forefoot pressure

Item #	Size	Quantity
6212	Left, One Size	1/pkg
6213	Right, One Size	1/pkg



Gel Bunion Shield

- Protects and helps relieve pain caused by bunions
- Gel dot cushions, protects and reduces pressure
- Mineral oil gel soothes and moisturizes the skin

Item #	Size	Quantity
6209	One Size	1/pkg





Gel Forefoot Pad

- Eases burning sensation in the ball of foot while preventing callus build-up
- Gel absorbs shock and gently soothes burning sensation in the metatarsal
- Built-in toe spreader helps reduce discomfort between toes while securing product in place

Item #	Size	Quantity
6211	One Size	2/pkg



Gel Ball of Foot Pad

- Absorbs shock and pressure on the ball of foot
- Gently massages the bottom of the foot
- Helps prevent calluses by reducing friction and abrasion

Item #	Size	Quantity
6214	One Size	2/pkg



Cushion Ball of Foot Pad

- Absorbs shock and pressure on the ball of foot
- Gently massages the bottom of the foot
- Foam provides lightweight comfort

Item #	Size	Quantity
6215	One Size	2/pkg





Heel Cups

- Relieves heel pain caused by plantar fasciitis and heel spurs
- Helps reduce lower back pain and discomfort
- Silicone construction provides all-day shock absorption

Item #	Size	Quantity
6210	One Size	2/pkg



Cushion Heel Pads

- Helps alleviate heel pain
- Self-adhesive backing for proper placement
- Reduces discomfort caused by plantar fasciitis

Item #	Size	Quantity
6221	Small	2/pkg
6222	Medium	2/pkg
6223	Large	2/pkg



Arch Support

- Provides gentle compression and comforting support to arch of foot
- Comfortably lifts and supports the plantar fascia
- Helps relieve heel and arch pain

Item #	Size	Quantity
6224	One Size	2/pkg





Prosthetics

Silipos gel products provide patients with an unmatched level of comfort and confidence. Unlike silicone, our mineral oil gels move with the skin and will not pull or irritate the residual limb. Moisturizing mineral oils are released from our gels which aids in the reduction of scar tissue and helps with overall skin conditioning. We offer a wide array of products that are sure to meet the needs of your practice and your patients.

Our prosthetic products are field tested with industry experts in order to assure patient comfort and satisfaction.

Silosheath™

Silosheath™ provides optimal range of motion, comfort, and the ability to be used with other interfaces as necessary. Made with Silipos mineral oil gel and nylon outer material, Silosheath™ is available in 3 different styles:



Silosheath™ Original

- 2ply, 2mm gel
- Packaged in a box

Item #	Size
13315	1
13325	2
13335	3
13345	4

Item #	Size
13355	5
13365	6
13375	7
13385	8



Silosheath™ Double Cushion

- 2 ply, 3mm gel, Extra gel at distal end
- Packaged in a box

Item #	Size
13415	1
13425	2
13435	3
13445	4

Item #	Size
13455	5
13465	6
13475	7
13485	8





Silosheath™ Extra Life

- 3 ply, 2mm gel, increased durability
- Packaged in a box

Item #	Size
13815	1
13825	2
13835	3
13845	4

Item #	Size
13855	5
13865	6
13875	7
13885	8

Silosheath™ Sizing Guidelines

Size	Description	Gel Length	Length	Distal Circumference	Proximal Circumference
					
1	Child	8" (20cm)	10"-13" (25cm-33cm)	6.5"-9" (17cm-23cm)	6.5"-10" (17cm-25cm)
2	Short / Narrow	10" (25cm)	15"-17" (38cm-43cm)	7"-10.5" (18cm-27cm)	7"-12" (18cm-30.5cm)
3	Short / Regular	10" (25cm)	15"-17" (38cm-43cm)	8"-12.5" (20cm-32cm)	8"-14" (20cm-35.5cm)
4	Medium / Narrow	13" (33cm)	17"-19" (43cm-48cm)	8"-12.5" (20cm-32cm)	8"-14" (20cm-35.5cm)
5	Medium / Regular	10" (25cm)	17"-19" (43cm-48cm)	9"-14" (23cm-35.5cm)	9"-16" (23cm-40.5cm)
6	Large / Narrow	13" (33cm)	18"-21" (46cm-53cm)	9"-14" (23cm-35.5cm)	9"-16" (23cm-40.5cm)
7	Large / Regular	13" (33cm)	18"-21" (46cm-53cm)	11"-16" (28cm-41cm)	11"-18" (28cm-46cm)
8	Large / Wide	13" (33cm)	18"-21" (46cm-53cm)	12"-19" (30.5cm-48cm)	12"-21" (30.5cm-53cm)

All sizes are approximate.

All Silosheath™ products offer the following benefits:

- Durability
- Conformability
- Maximum Protection
- Unmatched patient comfort and satisfaction

ComfortZone™ Gel Liners

Silipos Gel Liners offer the highest level in quality, protection and comfort. Whether your patient requires maximum protection or the added benefit of Silipos mineral oil gel, we have the product to fit your need. ComfortZone™ Gel Liners are available in 3 different styles:



Single Socket Gel Liner

- 3 ply, 3mm of exposed gel, CoolMax® exterior fabric, tapered shape
- Packaged in a box

Item #	Size	Ply
12725	Short / Narrow	3
12755	Short / Regular	3
12765	Medium / Regular	3



Double Socket Gel Liner

- 5 ply, 3mm gel, two layers of CoolMax® exterior fabric, tapered shape
- Packaged in a box

Item #	Size	Ply
12625	Short / Narrow	5
12655	Short / Regular	5
12665	Medium / Regular	5







Ultra Cushion Gel Liner

- 3 ply or 5 ply, 3mm gel, two layers of cotton material, nylon-reinforced distal end
- Packaged in a box

Item #	Size	Ply
12555	I - Short / Regular	3
12565	II - Medium / Regular	3
12575	III - Medium / Wide	3
12515	I - Short / Regular	5
12525	II - Medium / Regular	5
12535	III - Medium / Wide	5

CoolMax® is a registered trademark of Invista Corp.

ComfortZone™ Single Socket Gel Liner Sizing Guidelines

Item #	Size	Ply	Gel Length	Length of Liner	Toe Circumference	Top Circumference
						
12725	Short / Narrow	3	10" (25cm)	12" (30.5cm)	9"-12" (23cm-30.5cm)	15"-18" (38cm-45cm)
12755	Short / Regular	3	10" (25cm)	12" (30.5cm)	12"-14" (30.5cm-35.5cm)	15"-18" (38cm-45cm)
12765	Medium / Regular	3	12" (30.5cm)	16" (40.5cm)	12"-14" (30.5cm-35.5cm)	15"-18" (38cm-45cm)

All sizes are approximate.

ComfortZone™ Double Socket Gel Liner Sizing Guidelines

Item #	Size	Ply	Gel Length	Length of Liner	Toe Circumference	Top Circumference
12625	Short / Narrow	5	10" (25cm)	12" (30.5cm)	9"-12" (23cm-30.5cm)	15"-18" (38cm-45cm)
12655	Short / Regular	5	10" (25cm)	12" (30.5cm)	12"-14" (30.5cm-35.5cm)	15"-18" (38cm-45cm)
12665	Medium / Regular	5	12" (30.5cm)	16" (40.5cm)	12"-14" (30.5cm-35.5cm)	15"-18" (38cm-45cm)

All sizes are approximate.

ComfortZone™ Ultra Cushion Gel Liner Sizing Guidelines

Item #	Size	Ply	Gel Length	Length of Liner	Distal Circumference	Proximal Circumference
12555	I - Short / Regular	3	10" (25cm)	12.5" (32cm)	7"-11" (18cm-28cm)	8"-13" (20cm-33cm)
12515	I - Short / Regular	5	10" (25cm)	12.5" (32cm)	7"-11" (18cm-28cm)	8"-13" (20cm-33cm)
12565	II - Medium / Regular	3	12" (30.5cm)	15"-17" (38cm-43cm)	8"-12.5" (20cm-32cm)	8"-14" (20cm-35.5cm)
12525	II - Medium / Regular	5	12" (30.5cm)	17"-19" (43cm-48cm)	8"-12.5" (20cm-32cm)	8"-14" (20cm-35.5cm)
12575	III - Medium / Wide	3	12" (30.5cm)	17"-19" (43cm-48cm)	9"-14" (23cm-35.5cm)	9"-16" (23cm-40.5cm)
12535	III - Medium / Wide	5	12" (30.5cm)	18"-21" (46cm-53cm)	9"-14" (23cm-35.5cm)	9"-16" (23cm-40.5cm)

All sizes are approximate.

All ComfortZone™ Gel Liner products offer the following benefits:

- Prevent abrasion and skin breakdown
- Optimum cushion and protection
- Ease of mobility
- Unmatched patient comfort and satisfaction



SiloLiner™ is the most skin-friendly and easy to use solution for geriatric-transfer and post-op patients.

- Excellent for LTC patients
- Easy to Donn and Doff: Reinforced cuffs prevent liner from rolling down and irritating the skin
- Made with our softest gel that moisturizes and conditions skin grafts and scars to aid healing
- Outer material provides wide range of motion
- 6mm gel thickness
- Packaged in a box



SiloLiner™ Pin Liner, Sizing Guidelines

Item #		Flat Width	Distal Circumference at urethane attachment	Proximal Circumference MPT
1460	1464	3\" (7.6cm)	8\"-12\" (20cm-30.5cm)	8\"-14\" (20cm-35.5cm)
1461	1465	4\" (10.2cm)	9.5\"-13.5\" (24cm-34cm)	9.5\"-16\" (24cm-40.5cm)
1462	1466	5\" (12.7cm)	12\"-15\" (30.5cm-38cm)	12\"-20\" (30.5cm-51cm)
1463	1467	6\" (15.2cm)	14\"-16\" (35.5cm-40.5cm)	14\"-26\" (35.5cm-66cm)

All sizes are approximate.

Pin Liner

Item #	Length	Flat Width
1460	12" (30.5cm)	3" (7.6cm)
1461	12" (30.5cm)	4" (10.2cm)
1462	12" (30.5cm)	5" (12.7cm)
1463	12" (30.5cm)	6" (15.2cm)
1464	14" (35.5cm)	3" (7.6cm)
1465	14" (35.5cm)	4" (10.2cm)
1466	14" (35.5cm)	5" (12.7cm)
1467	14" (35.5cm)	6" (15.2cm)



SiloLiner™ Cushion Liner, Sizing Guidelines

Item #		Flat Width	Distal Circumference at urethane attachment	Proximal Circumference MPT
1612	1622	3" (7.6cm)	8"-14" (20cm-35.5cm)	8"-14" (20cm-35.5cm)
1613	1612	4" (10.2cm)	9.5"-16" (24cm-40.5cm)	9.5"-16" (24cm-40.5cm)
1614	1624	5" (12.7cm)	12"-17" (30.5cm-43cm)	12"-20" (30.5cm-51cm)
1615	1625	6" (15.2cm)	14"-19" (35.5cm-48cm)	14"-26" (35.5cm-66cm)

All sizes are approximate.

Cushion Liner

Item #	Length	Flat Width
1612	12" (30.5cm)	3" (7.6cm)
1613	12" (30.5cm)	4" (10.2cm)
1614	12" (30.5cm)	5" (12.7cm)
1615	12" (30.5cm)	6" (15.2cm)
1622	14" (35.5cm)	3" (7.6cm)
1623	14" (35.5cm)	4" (10.2cm)
1624	14" (35.5cm)	5" (12.7cm)
1625	14" (35.5cm)	6" (15.2cm)

Low Activity Liner

The Low Activity Liner provides superior skin conditioning and comfort for low to moderate activity level patients.

- Lined with a medium durometer gel
- Stabilizing Matrix reduces distal elongation and minimizes pistoning
- Liner easily conforms to limbs
- Excellent for sensitive skin
- Retrofittable to Silipos Explorer™ sockets
- Available in 3mm, 6mm, and 9mm gel thickness
- Packaged in a box



Pin Liner, 3mm Gel Thickness

Item #	Size	Style
18101	Medium	Uniform
18107	Medium +	Uniform
18113	Large	Uniform
18116	Large +	Uniform
18119	XLarge	Uniform

Item #	Size	Style
18102	Medium	Tapered
18108	Medium +	Tapered

Pin Liner, 6mm Gel Thickness

Item #	Size	Style
18103	Medium	Uniform
18109	Medium +	Uniform
18114	Large	Uniform
18117	Large +	Uniform
18120	XLarge	Uniform



Item #	Size	Style
18104	Medium	Tapered
18110	Medium +	Tapered

Pin Liner, 9mm Gel Thickness

Item #	Size	Style
18105	Medium	Uniform
18111	Medium +	Uniform
18115	Large	Uniform
18118	Large +	Uniform

Item #	Size	Style
18106	Medium	Tapered
18112	Medium +	Tapered

Low Activity Liner Sizing Guidelines

	Proximal Circumference		Distal Circumference	
	 Measure 15" from Distal End		 Measure 1.5" from Distal End	
Size	Relaxed	Fully Stretched	Relaxed	Fully Stretched
Medium	9" (23cm)	16" (41cm)	8" (20cm)	12" (30cm)
Medium+	13" (33cm)	21" (53cm)	9" (23cm)	12" (30cm)
Large	15" (38cm)	22" (56cm)	9" (23cm)	15" (38cm)
Large+	16" (41cm)	27" (68cm)	12" (30cm)	16" (41cm)
X-Large	16" (41cm)	30" (76cm)	12" (30cm)	20" (51cm)

All sizes are approximate.



Cushion Liner, 3mm Gel Thickness

Item #	Size	Style
18141	Medium	Uniform
18147	Medium +	Uniform
18153	Large	Uniform
18156	Large +	Uniform
18159	XLarge	Uniform

Item #	Size	Style
18142	Medium	Tapered
18148	Medium +	Tapered

Cushion Liner, 6mm Gel Thickness

Item #	Size	Style
18143	Medium	Uniform
18149	Medium +	Uniform
18154	Large	Uniform
18157	Large +	Uniform
18160	XLarge	Uniform

Item #	Size	Style
18144	Medium	Tapered
18150	Medium +	Tapered

Cushion Liner, 9mm Gel Thickness

Item #	Size	Style
18145	Medium	Uniform
18151	Medium +	Uniform
18155	Large	Uniform
18158	Large +	Uniform

Item #	Size	Style
18146	Medium	Tapered
18152	Medium +	Tapered



DuraGel™ offers an unparalleled combination of durability, comfort, and skin protection for low to high activity level patients.

- Made with a durable gel formulation designed to help increase liner life and maintain comfort
- Ideal for moderately active patients engaging in every day living activities
- Nylon outer fabric engineered to provide superior abrasion resistance and gel bond
- Available in common pin and cushion profiles: 3mm, 6mm, and 9mm; Tapered and Uniform styles
- Packaged in a box



Pin Liner, 3mm Gel Thickness

Item #	Size	Style
18211	Medium	Uniform
18212	Medium +	Uniform
18213	Large	Uniform
18214	Large +	Uniform
18215	X-Large	Uniform

Item #	Size	Style
18216	Medium	Tapered
18217	Medium +	Tapered
18218	Large	Tapered

Pin Liner, 6mm Gel Thickness



Item #	Size	Style
18201	Medium	Uniform
18202	Medium +	Uniform
18203	Large	Uniform
18204	Large +	Uniform
18205	X-Large	Uniform

Item #	Size	Style
18206	Medium	Tapered
18207	Medium +	Tapered
18208	Large	Tapered

Pin Liner, 9mm Gel Thickness

Item #	Size	Style
18221	Medium	Uniform
18222	Medium +	Uniform
18223	Large	Uniform
18224	Large +	Uniform

DuraGel™ Sizing Guidelines

Proximal Circumference			Distal Circumference	
 Measure 15" from Distal End			 Measure 1.5" from Distal End	
Size	Relaxed	Fully Stretched	Relaxed	Fully Stretched
Medium	9" (23cm)	15" (38cm)	8" (20cm)	13" (33cm)
Medium+	13" (33cm)	20" (51cm)	9" (23cm)	13" (33cm)
Large	15" (38cm)	21" (53cm)	9" (23cm)	16" (41cm)
Large+	16" (41cm)	26" (66cm)	12" (30cm)	17" (43cm)
X-Large	16" (41cm)	29" (74cm)	12" (30cm)	21" (53cm)

All sizes are approximate.



Cushion Liner, 3mm Gel Thickness

Item #	Size	Style
18251	Medium	Uniform
18252	Medium +	Uniform
18253	Large	Uniform
18254	Large +	Uniform
18255	X-Large	Uniform

Item #	Size	Style
18256	Medium	Tapered
18257	Medium +	Tapered

Cushion Liner, 6mm Gel Thickness

Item #	Size	Style
18241	Medium	Uniform
18242	Medium +	Uniform
18243	Large	Uniform
18244	Large +	Uniform
18245	X-Large	Uniform

Item #	Size	Style
18246	Medium	Tapered
18247	Medium +	Tapered
18248	Large	Tapered

Cushion Liner, 9mm Gel Thickness

Item #	Size	Style
18261	Medium	Uniform
18262	Medium +	Uniform
18263	Large	Uniform
18264	Large +	Uniform



DuraGel™ Suspension Sleeve

The DuraGel™ Suspension Sleeve combines the same gel used in our DuraGel™ liners with a flexible outer fabric to provide improved durability and range of motion.

- 3mm gel thickness and outer fabric offers easy knee flexion
- Lined with durable gel formulation for increased durability
- Mineral oil gel conditions and moisturizes the skin
- Packaged in a polybag



Item #	Size
18181	Medium
18182	Medium +
18185	Medium+ Distal+

Item #	Size
18183	Large
18184	Large +



DuraGel™ Reinforced Suspension Sleeve

Internal fabric reinforcement provides increased durability without limiting flexibility. This sleeve maintains suspension and provides added security for its wearer.

- Packaged in a polybag



Item #	Size
18191	Medium
18192	Medium +
18195	Medium+ Distal+

Item #	Size
18193	Large
18194	Large +

DuraGel™ Suspension Sleeve Sizing Guidelines

Size	Proximal Circumference		Distal Circumference	
	Relaxed	Fully Stretched	Relaxed	Fully Stretched
Medium	9" (23cm)	15" (38cm)	8" (20cm)	13" (33cm)
Medium+	14.25" (36cm)	26.75" (68cm)	9.25" (23cm)	16" (40cm)
Medium+ Distal+	14.25" (36cm)	26.75" (68cm)	9.75" (24.5cm)	17" (43cm)
Large	13.5" (34cm)	24" (61cm)	12" (31cm)	18" (46cm)
Large+	16.25" (41cm)	28" (71cm)	11.25" (29cm)	20" (51cm)

All sizes are approximate.



BK Suspension Sleeve

- Provides maximum suspension
- Lightweight and durable
- Ease of application
- Seamless construction
- Made with our softest gel that moisturizes and conditions the skin
- Packaged in a polybag



HD-BK Suspension Sleeve

- Helps to prevent rotation
- Internal fabric reinforcement increases durability
- Will not bunch-up behind the popliteal area
- Ease of application
- Made with our softest gel that moisturizes and conditions the skin
- Packaged in a polybag



Arm Gel Suspension Sleeve

- Compatible with myo-electric and cosmetic prosthesis
- May also be used as BK sleeve for children and small adults
- Lined with a soft mineral oil gel moisturizes and conditions the skin
- Packaged in a polybag



Item #	Size	Approx. Length	Approx. Circumference
13055	I	12" (30.5cm)	11"-17" (28cm-43cm)
13065	II	14" (35.5cm)	13"-20" (32.5cm-51cm)
13075	III	16" (40.5cm)	14"-22" (35.5cm-56cm)
13085	IV	16" (40.5cm)	15"-28" (38cm-71cm)

Item #	Size	Approx. Length	Approx. Circumference
1310	I	12" (30.5cm)	11"-17" (28cm-43cm)
1311	II	14" (35.5cm)	13"-20" (32.5cm-51cm)
1312	III	16" (40.5cm)	14"-22" (35.5cm-56cm)





Item #	Approx. Length	Approx. Circumference
13005	8-13" (20-30cm)	8.75" (22cm)



Symes Gel Liner

- 5 ply, 3mm gel
- Retrofittable to most Syme's socket designs
- Prevents abrasion and skin breakdown
- Optimum conformability
- Made with mineral oil gel
- CoolMax® Material provides elasticity, comfort, and breathability
- Packaged in a polybag

Item #	Size
1302	One Size Fits Most

Item #	Length	Gel Length	Distal Circumference	Proximal Circumference
				
1302	16.5" (42cm)	15" (38cm)	10"-14.5" (25cm-37cm)	15"-20" (38cm-51cm)

All sizes are approximate.



Partial Foot Sock

- For use with transmetatarsal amputees
- CoolMax® Material provides elasticity, comfort, and breathability
- May be used with foot orthoses or shoes with fillers
- Available in 2 sizes
- Mineral oil gel covers site of amputation
- Packaged in a polybag

Item #	Size
10805	Regular
10815	Extra Long

CoolMax® is a registered trademark of Invista Corp.



Mineral Oil Gel Roll-On

- Apply to distal end of limb for added comfort and protection
- Helps to enhance socket liner fit
- Adheres to Silipos gel liners to ease application
- Can be cut to size
- Fully constructed of mineral oil gel to moisturize and protect skin

Item #	Size	Gel Thickness	Approx. Outer Circ.	Approx. Gel Length
1495	Narrow	2mm	7.25" (18.5cm)	6" (15cm)
1496	Medium	2mm	8.5" (21.5cm)	6" (15cm)
1497	Wide	2mm	9.25" (23.5cm)	6" (15cm)



Distal End Pad

- Provides extra protection for shock and abrasion at the distal end
- Aids in enhancing socket fit for limbs undergoing changes in volume or mass
- Can be cut to size
- Packaged in a polybag

Item #	Size
12005	One Size



Skin Care

Silipos gels cushion and protect specific parts of the body against impact, pressure, and shear forces, while simultaneously releasing moisture-rich mineral oil that gently hydrates, softens, and conditions the skin without the use of additional creams or lotions.

Visco-elastic properties enable our gels to move synergistically with the skin which prevents abrasion and discomfort. Our user-friendly products easily adapt to a wide range of body types and builds to ensure high patient compliance.

All of our gels are clear, non-toxic, and hypoallergenic. Technologically advanced manufacturing capabilities allow Silipos to create a vast assortment of gels for a variety of uses. Additives such as vitamins and minerals, antioxidants, moisturizers, and fragrances can be applied to our gels to address specific concerns across numerous markets. Silipos gels are washable and reusable, and with proper care will continue to release skin-beneficial mineral oils throughout the life of the product.



Gel Therapy Gloves & Socks

- Mineral oil gel comfortably conforms to the hands or feet to continuously moisturize and condition the skin
- Beneficial for skin conditions such as eczema, psoriasis, and dermatitis
- Designed with non-restrictive cuff
- Can be used as warming gloves/socks for arthritic conditions
- Gel is hypoallergenic and does not contain any additives, fragrances, or latex
- Gloves and socks can be worn at night for intensive treatment (products are non-ambulatory)

Gel Therapy Gloves *(Shrinkwrap packaged with backerboard)*

Item #	Size	Color	Quantity
786	One Size	<input type="checkbox"/> White	1 pair/pkg

Gel Therapy Socks *(Packaged in a box)*

Item #	Size	Color	Quantity
788	One Size	<input type="checkbox"/> Oatmeal	1 pair/pkg



NouveauDerm™ Gel Gloves & Socks

- Mineral oil gel comfortably conforms to the hands or feet to continuously moisturize and condition the skin
- Use as treatment for dry, hard, or callused skin
- Infused with an assortment of fragrances and available in 3 different colors
- Packaged in a box

NouveauDerm™ Gel Gloves

Item #	Size	Color	Scent	Quantity
206	One Size	<input type="checkbox"/> Yellow	Green Tea	1 pair/pkg
207	One Size	<input type="checkbox"/> Green	Green Tea	1 pair/pkg
208	One Size	<input type="checkbox"/> Oatmeal	Green Tea	1 pair/pkg

NouveauDerm™ Gel Socks

Item #	Size	Color	Scent	Quantity
216	One Size	<input type="checkbox"/> Yellow	Green Tea	1 pair/pkg
217	One Size	<input type="checkbox"/> Green	Green Tea	1 pair/pkg
218	One Size	<input type="checkbox"/> Oatmeal	Green Tea	1 pair/pkg



Gel-Care™


- 2mm thick layer of self-adhesive mineral oil gel with conforming material backing
- Moisturizes the skin and helps to reduce scar tissue
- Hypoallergenic, washable and reusable
- Cost effective solution for scar management
- Packaged in an envelope-style box

Item #	Size	Quantity
760	4" x 4" (10cm x 10cm)	1/pkg
761	5.5" x 8" (14cm x 20cm)	1/pkg
762	1.5" x 10" (4cm x 25cm)	2/pkg



Post Surgical Breast Forms

- 0.5mm thick layer of self-adhesive mineral oil gel with conforming nylon backing
- For use after breast augmentation procedures
- Allows patient to wear clothing next to hyper sensitive skin areas
- Packaged in a reusable tray




Item #	Size	Quantity
740	Small, 32-36 (81.2cm x 91.4cm)	1 pair/pkg
741	Large, 36-42 (91.4cm x 106.7cm)	1 pair/pkg



Post Surgical Abdominal Form

- 0.5mm thick layer of self-adhesive mineral oil gel with conforming nylon backing
- For use after abdominoplasty procedures
- Allows patient to wear clothing next to hyper sensitive skin areas
- Packaged in a reusable tray



Item #	Size	Quantity
744	Small, 17.3" x 2.4" (44cm x 6cm)	1 pair/pkg
745	Large, 17.7" x 3.5" (45cm x 9cm)	1 pair/pkg



THE GLOBAL LEADER IN GEL TECHNOLOGIES®

Silipos, Inc.
7049 Williams Road
Niagara Falls, NY 14304
United States
Toll Free: 800-229-4404
Phone: +1 (716) 283-0700
Fax: +1 (716) 283-0600

Langer UK
Brookhouse Way
Cheadle
Stoke-on-Trent
Staffordshire, England, ST10 1SR
Phone: +44 (0) 845-658-1650
Fax: +44 (0) 845-658-1652

Silipos Asia Pacific
P.O. Box 5122
Mordialloc
Victoria 3195
Australia
Phone: +61 (3) 9590-9818
Fax: +61 (3) 9590-9003

www.silipos.com